

Product Awareness Bulletin – April 2018

PRODUCT UPDATES: Pride Jazzy Select Elite 2 & Go Chair NG

Pride Go Chair NG (Next Generation)

Key features:

- The re-engineered Go Chair offers enhanced performance and comfort, along with feather-touch disassembly. With an increased weight capacity of 136 kg and a sleek, bold look, the new Go Chair makes travel easy
- Portable travel design – the folding seat and lightweight frame easily disassemble into five pieces for convenient storage or transportation (Designed to fit into the trunk of most vehicles)
- Small turning radius and compact design - provide superior handling ideal for confined spaces
- 5-point swivel seat with latch release - facilitates user transfers
- Height & width adjustable armrests
- Dual swing-away, under seat storage bins
- 60" lap belt
- Easy-to-remove battery pack – offers the convenience of indoor, off-board charging
- Colours available: Red
- Target user - Users with low/basic seating and positioning needs looking for an economical power wheelchair for indoor use

	
<i>Product Code:</i>	<i>PWS654200RE</i>
<i>Maximum User Weight:</i>	<i>136kg</i>
<i>Base Weight:</i>	<i>30kg (front & rear)</i>
<i>Seat Weight:</i>	<i>13.6kg</i>
<i>Battery Weight:</i>	<i>14kg (battery pack)</i>
<i>Maximum Speed:</i>	<i>Up to 5.95km/h</i>
<i>Maximum Range Per Charge:</i>	<i>Up to 14km</i>
<i>Turning Radius:</i>	<i>64.5cm (25.4") with foot platform</i>
<i>Standard Seat Dimensions:</i>	<i>46 x 43cm (18 x 17")</i>
<i>Length/Width:</i>	<i>84cm (33")/57cm (22.6")</i>
<i>Suspension Type:</i>	<i>Limited*</i>

Click [here](#) to view the product online

Pride Jazzy Select Elite 2

Key features:

- **Exclusive to Aidacare**
- In-line, front wheel drive technology for stability and manoeuvrability both indoors and outdoors
- Easy access to batteries and simple main frame design for ease of serviceability
- Comfort, high-back seat with headrest
- Off board charging
- Colours available: Jazzy Red
- Target User - Users who want a chair for indoor and outdoor use that is a blend of power, performance and style

	
<i>Product Code:</i>	<i>PWS654255</i>
<i>Maximum User Weight:</i>	<i>136kg</i>
<i>Base Weight:</i>	<i>51kg</i>
<i>Seat Weight:</i>	<i>17.7kg (comfort high back)</i>
<i>Battery Weight:</i>	<i>11.1kg each</i>
<i>Maximum Speed:</i>	<i>Up to 6.44km/h</i>
<i>Maximum Range Per Charge:</i>	<i>Up to 22km</i>
<i>Turning Radius:</i>	<i>61cm (24")</i>
<i>Standard Seat Dimensions:</i>	<i>50.8 x 50.8cm (20 x 20")</i>
<i>Length/Width:</i>	<i>87.3cm (37") x 58.4cm (23")</i>
<i>Suspension Type:</i>	<i>Limited</i>

Click [here](#) to view the product online