





# aspire LIFT RECLINE CHAIRS RANGE OVERVIEW


	GOLD		SILVER			BRONZE				
	RAPHAEL QUATTRO LIFT RECLINE CHAIR	DA VINCI QUATTRO LIFT RECLINE CHAIR	SIGNATURE SPACE SAVER LIFT RECLINE CHAIR	SIGNATURE DUAL ACTION LIFT RECLINE CHAIR		PORTO SPACE SAVER PETITE LIFT RECLINE CHAIR	VERONA SPACE SAVER LIFT RECLINE CHAIR	MAUI DUAL ACTION LIFT RECLINE CHAIR	MONTANA SINGLE ACTION MAXI LIFT RECLINE CHAIR	
										
Description	4 Motor - Quattro	4 Motor - Quattro	1 Motor - Space Saver	2 Motor - Dual Action		1 Motor - Space Saver	1 Motor - Space Saver	2 Motor - Dual Action	1 Motor - Heavy Duty	
Functions										
Size	Medium	Medium	Small	Medium	Large	Petite	Medium	Medium	Maxi	
Upholstery	Microsuede Leather-look	Microsuede	Warwick Global Suede & Vinyl			Soft Suede	Quilted Upholstery	Quilted Upholstery	Microsuede Leather-look	
Seat to Floor	480 mm	480 mm	470 mm	480 mm	510 mm	470 mm	480 mm	480 mm	510 mm	
Seat Depth	510 mm	510 mm	450 - 475 mm	490 - 515 mm	510 - 535 mm	500 mm	500 mm	500 mm	530 mm	
Seat Width	530 mm	530 mm	450 mm	480 mm	500 mm	490 mm	560 mm	560 mm	630 mm	
Overall Length	950 mm	950 mm	1010 mm	1050 mm	1110 mm	900 mm	800 mm	800 mm	880 mm	
Overall Width	810 mm	810 mm	780 mm	780 mm	780 mm	770 mm	820 mm	820 mm	940 mm	
Backrest Height	620 mm	630 mm	590 mm	620 mm	670 mm	700 mm	630 mm	630 mm	660 mm	
SWL	130 kg	159 kg	145 kg	155 kg		150 kg	130 kg	130 kg	220 kg	
AUST R ARTG No.	306075	216484	230186	230186		306075	306075	306075	306075	

## COLOUR RANGE



### ASPIRE ACCESSORIES

**Chair Protector**  
Universal Size  
CHP227600



**Signature Chairs - Head & Armrest Cover Set**

