

# NSW HEALTH STE830 EQUIPMENT CATALOGUE

**NOVEMBER 2019** 

VERSION: 1 UPDATED: NOV 2019



Quality Equipment Enterprise Solutions Home & Community





# Aidacare & NSW Health

Aidacare are proud to be recognised as a NSW HealthShare / EnableNSW contracted service provider for a large range of products included in the recently announced Beds and Mattresses contract. Our commitment to all NSW HealthShare stakeholders including prescribing health care professionals and most importantly, the NSW HealthShare client, is to provide the very best service and product quality available.

For the healthcare prescriber, Aidacare offer industry leading support to assist in the provision of equipment for your clients.

- All contracted products are available for trial, anywhere in New South Wales, to enable assessment of suitability for your clients specific needs. Our equipment consultants can meet with you at any location, including client's private homes, to demonstrate and trial items.
- Dedicated NSW HealthShare display stock located in each of the Aidacare showrooms is available to enable prescribers to meet with clients on-site to trial and assess all required items in a professional environment.
- Many Aidacare showrooms even offer prescribers use of a private office and meeting room to utilise while meeting with clients or just to catch up on admin and reports.
- Specialist Complex Rehab & Seating staff to assist and support with prescription of customised and specialised equipment.
- Clinical education and support via the Aidacare Professional Training program and workshops.
- Service and repair of all contracted items with our workshop and mobile technicians.

General customer service and support including the booking of trials and any general advice and information required.

Our electronic catalogue is available for download at www.aidacare.com.au/globalassets/resources

For information on any of the extensive range of products available from Aidacare, please browse the complete online Aidacare catalogue at www.aidacare.com.au

For additional information, please contact your local Aidacare office directly on 1300 133 120 or email nswhealth@aidacare.com.au

EnableNSW provides equipment and services to people in NSW with chronic health conditions or disability to assist them with mobility, communication and self-care. We also provide assistance for people who need to travel significant distances for specialist medical care.

Aids and Equipment Program - The Aids and Equipment Program is a NSW Government initiative to assist eligible residents of NSW, who have a life-long or long-term disability, to live and participate within their community by providing appropriate equipment, aids and appliances.

Information on eligibility for the program, the range of equipment available, the lodgement of applications, mutual responsibilities and expectations, and steps for complaints and appeals is outlined in the Policy Directive - "Assistive Technology for Communication, Mobility, Respiratory Function & Self-Care". https://www1.health.nsw.gov.au/PDS/pages/doc.aspx?dn=PD2011\_027

People with a disability interested in seeking further information about the program should discuss their query with a healthcare provider, or Enable NSW on 1800 ENABLE (1800 362 253)

As part of HealthShare NSW, EnableNSW is responsible for the administration of NSW Health disability support and other assistance programs including:

- Aids and Equipment Program
- Continence Assistance
- Home Respiratory Program (HRP)
- Prosthetic Limb Service (formerly known as ALS)
- Specialised Equipment Essential for Discharge (SEED)

For more information, visit the EnableNSW website: http://www.enable.health.nsw.gov.au/

# **Booking a Service or Repair**

Our mobile-repair technicians can come to you to repair or maintain your equipment. Our qualified staff are industry trained, test and tag certified and Police checked.

#### How to request a repair:

Speak to your consultant to arrange your service & repairs.

## **General Equipment Enquiries:**

General Contact Details: 1300 133 120

NSW Health Email Enquiries: nswhealth@aidacare.com.au

#### **Service Department Regular Business Hours:**

Phone: 1300 133 120

Email: nswhealth@aidacare.com.au

Hours: 9:00am - 5:00pm

# **Drying Racks**




Stainless Steel 4 Bedpan 4 Urine Bottle Drying Rack UPA200270



Stainless Steel 6 Bedpan 6 Urine Bottle Drying Rack UPA200290

# **Utensil Trays**



Optima 3 / Maxi UPA200250



Optima 2 / Midi UPA200240

# Sanitisers & Macerators









Easier Loading

cycle start

opening

## Panamatic Optima 2

The Panamatic Optima 2 Washer Sanitiser (hot & cold water version of the Midi) has raised the standard to help meet the ongoing hygiene challenges in Hospital & Aged Care environments.

#### Features:

- Improved wash performance
- Efficient water and power use
- Hands-free open, close and start
- MicrobeSafe antimicrobial lid
- Steam disinfection cycle

#### **Optional Accessories**

UPA200270 - Bed Pan Rack w/ Drip Tray (4 Pan / 4 Bottle)

UPA200290 -Bed Pan Rack w/ DripTray (6 Pan / 6 Bottle)

UPA200240 - Utensil tray + sides (2 piece set)

Code	UPM100200
Electrical supply	Single Phase 230 V, 13 AMP
Electrical output/input	Motor - 0.75 KW, Heater - 2.5KW
Dimensions (w x d x h)	500 x 600 x 1042 mm
Water supply	15 mm Hot & Cold 1 - 6 bar
Soil connection	110 mm connections supplied
Operation cycle time	Approx 5 mins / 30ml EcoWash
Steam disinfection	Exceeds Ao 700, 85° x 75 secs
Typical capacity	2 bedpans + 2 urine bottles/cycle
Litres	63.4 Litres
Warranty	1 Year
Description	Clinical Monitoring System
Code	UPA200300



Easier Loading Greater Versatility

Hands free

## Panamatic Optima 3

The Panamatic Optima 3 Washer Sanitiser is the hot & cold water version of the Maxi. The large wash chamber accommodates a range of utensils and urine bottles per cycle.

#### Features:

- Improved wash performance
- Hands-free open, close and start
- Large capacity chamber

#### **Optional Accessories**

UPA200270 - Bed Pan Rack w/ DripTray (4 Pan / 4 Bottle)

UPA200290 - Bed Pan Rack w/ DripTray (6 Pan / 6 Bottle)

UPA200250 - Utensil tray + sides (2 piece set)

Code	UPM100250
Electrical supply	Single Phase 230 V, 13 AMP
Electrical output/input	Motor - 0.75 KW, Heater - 2.5 KW
Dimensions (w x d x h)	600 x 600 x 900 mm
Water supply	15 mm Hot & Cold water 1 - 6 bar
Soil connection	110 mm connections supplied
Operation cycle time	Approx 5 mins - 30ml Ecowash
Steam disinfection	Exceeds Ao 700, 85° x 75 secs
Typical capacity	3 bedpans + 3 urine bottles/cycle
Litres	71.9 Litres
Warranty	1 Year
Description	Clinical Monitoring System
Code	UPA200310



Easier Loading Greater Versatility

opening

## Panamatic Maxi

The Panamatic Maxi Washer Sanitiser is the cold water version of the Optima 3. The large wash chamber accommodates a range of utensils and urine bottles per cycle.

#### Features:

- Cold water wash
- Large chamber capacity
- Steam disinfection cycle

#### Optional Accessories

UPA200270 - Bed Pan Rack w/ DripTray (4 Pan / 4 Bottle)

UPA200290 - Bed Pan Rack w/ DripTray (6 Pan / 6 Bottle)

UPA200250 - Utensil tray + sides (2 piece set)

Code	UPM100350
Electrical supply	Single Phase 230 V, 13 AMP
Electrical output/input	Motor - 0.75 KW, Heater - 2.5 KW
Dimensions (w x d x h)	600 x 600 x 900 mm
Water supply	15 mm Cold water 0.5 - 10 bar
Soil connection	110 mm connections supplied
Operation cycle time	Approx 5 mins - 30ml Ecowash
Steam disinfection	Exceeds Ao 700, 85° x 75 secs
Typical capacity	3 bedpans + 3 urine bottles/cycle
Litres	71.9 Litres
Warranty	1 Year
Description	Clinical Manitoring System
Description	Clinical Monitoring System
Code	UPA200310



Easier Loading

Hands free cycle start

## Panamatic Midi

The Panamatic Midi Washer Sanitiser is the cold water version of the Optima 2.

#### Features:

- Cold water wash
- Efficient water and power use
- Hands-free open, close and start
- MicrobeSafe antimicrobial lid
- Large chamber capacity

#### **Optional Accessories**

UPA200270 - Bed Pan Rack w/ DripTray (4 Pan / 4 Bottle)

UPA200290 - Bed Pan Rack w/ DripTray (6 Pan / 6 Bottle)

UPA200240 - Utensil tray + sides (2 piece set)

Code	UPM100300
Electrical supply	Single Phase 230 V, 13 AMP
Electrical output/input	Motor - 0.425KW, Heater - 2.5KW
Dimensions (w x d x h)	500 x 600 x 1042 mm
Water supply	15 mm Cold water 0.5 - 10 bar
Soil connection	110 mm connections supplied
Operation cycle time	Approx 5 mins / 30ml EcoWash
Steam disinfection	Exceeds Ao 700, 85° x 75 secs
Typical capacity	2 bedpans + 2 urine bottles/cycle
Litres	63.4 Litres
Warranty	1 Year
Description	Clinical Monitoring System
Code	UPA200300





## THE DDC RANGE

# **CROSS-INFECTION CONTROL**

### **BREAKING THE CHAIN OF INFECTION · · · ·**



Infection Control



Water Efficient



Budget Friendly



Energy



Nationally Accredited

WASHER SANITISERS: STEAM DISINFECTION TO > AO 700

#### HOT & COLD WATER SYSTEMS



Panamatic
Optima 3













Hands Free Opening



Panamatic
Optima 2













Hands Free Open, Close & Cycle Start

#### COLD WATER SYSTEMS



Panamatic Maxi









Capacity:

3 Bed Pans &



Hands Free Opening



Panamatic Midi L









Hands Free Open, Close & Cycle Start

# State Wide Distribution Network

MOOREBANK
3A / 1 Moorebank Ave
Moorebank NSW 2170
Tel: (02) 8706 2300
Fax: (02) 9618 5111

SEVEN HILLS
Unit 2/4 Stanton Rd
Seven Hills NSW 2147
Tel: (02) 8706 2330
Fax: (02) 9624 4256

BLACKTOWN 1/34 Campbell St Blacktown NSW 2178 Tel: (02) 9621 4077

WOLLONGONG
140 Church St
Wollongong NSW 2500
Tel: (02) 8706 2330

ALBURY 347 Urana Rd Lavington NSW 2641 Tel: (02) 6040 3444 Fax: (02) 6040 3222

ORANGE
Unit 1/12-16 Ash St
Orange NSW 2800
Tel: 0437 134 295

TAMWORTH
Unit 3/248 Peel St
Tamworth NSW 2340
Tel: (02) 6761 2018
Fax: (02) 6761 2090

PORT MACQUARIE
171 Lake Rd
Port Macquarie, NSW, 2444
Tel: (02) 6581 0018
Fax: (02) 6581 0038

COFFS HARBOUR 2/177 Orlando St Coffs Harbour NSW 2450 Tel: (02) 6651 5333

BALLINA
18 Endeavour Cl
Ballina NSW 2478
Tel: (02) 6686 6644
Fax: (02) 6686 9383



CALL 1300 133 120

Call this number to be automatically diverted to your local Aidacare office (for the price of a local call).









