

TRANSPORT RAMP

Aspire for... Access, Mobility & Support



Durable all weather surfaces

FEATURES & BENEFITS

- Lightweight yet robust. A super lightweight composite fibre core is supported by an aluminium frame to support loads of up to 400kg.
- Anti-slip all weather surface.
- Raised tapered side walls keep wheels tracking along the ramp surface.
- Folding design with carry handle and fastening clip allows the ramp to be easily carried, manoeuvred into position or stowed when not in use.
- Striking orange colour provides high visibility and contrast.



Folding for quick access & storage



Caters to Bariatric Users with a Safe Working Load (SWL) of 400kg



2 YEAR WARRANTY

CODE	DESCRIPTION	LENGTH	TOTAL WIDTH	PLATFORM WIDTH	FOLDED WIDTH	FOLDED HEIGHT	PRODUCT WEIGHT
RAM656800	Aspire Transport Ramp	1500 mm	850 mm	800 mm	435 mm	55 mm	10 kg