

# A200S STANDING LIFTER

The A200S Standing Lifter optimises patient comfort through its depth-adjustable yoke system and lightweight aluminium construction. Designed for versatile home-to-facility use, this compact mobility solution features adjustable padded leg rests and intelligent diagnostics, enabling personalised standing assistance across diverse care environments.



#### Depth-adjustable Yoke

Yoke is depth adjustable for optimal patient fit.



**5** YEAR WARRANTY  
Frame

**3** YEAR WARRANTY  
Electronics

**1** YEAR WARRANTY  
Battery, Handset, and Cables

SWL  
**200kg**



### Intelligent Diagnostics

State of the art Linak Jumbo Care control system. Tracks usage and service requirements for safety and compliance.



### Build Quality

Strong powder-coated aluminium frame with industry leading Linak electronics, Tente castors and warranty provision.



### Depth-adjustable Yoke

Yoke is depth adjustable for optimal patient fit.




### Adjustable Padded Legrests

Provide comfort and support during transfers and can be adjusted to suit the patient.




### Ergonomic Handles

Provide multiple grip points to suit carers of different heights and strengths.



### Smart Handset

The latest slimline technology from Linak fits neatly in the hand with an intuitive control panel. Battery and service indicator lights assist with service scheduling.



### TENTE

2 x 75mm front castors for maximum under bed clearance and accessibility.

2 x 75mm rear-locking castors for safety and manoeuvrability.



### Safety


Emergency stop button, emergency lowering and anti-entrapment technology, with both patient and caregiver safety at heart.



### Powered Leg Function

Opens and closes the legs for lifting stability or accessibility.

## ACCESSORY



### Sling Compatibility

Compatible with a range of Aspire slings.



See our Sling Fitting and Safe Use Guide for full details

Code	Frame Material	Lifting Range	Under Bed Clearance	Turning Diameter	Base Length	External Base Width	Internal Base Width	Overall Height	Lifter Weight	SWL
LSS390640	Aluminium	692 - 1770	110	1160	1070	590 - 1040	519 - 986	1480	33kg	200kg

All measurements are in millimetres unless otherwise specified.