

aidacare INLINE WEIGH SCALE







BRANDS



PRODUCT CODE:

Year

Warrant v

LSS390680

Aidacare Aspire In-Line Weigh Scale

- Designed to be fitted on Aidacare Aspire Patient Lifters, however adaptable to most Lifter brands
- Allows for an accurate weight measurement to be taken by a carer when a patient is suspended in a sling
- Large weight capacity of up to 320kg
- Easy to read LCD display

SUITS AIDACARE ASPIRE PATIENT LIFTERS:



LSS390600 - A150F Aluminum Folding Lifter



LSS390610 - A205 Aluminium Patient Lifter



LSS390620 - 200 Steel Patient Lifter



LSS390630 - A320 Aluminium Bariatric Lifter

