




REHAB RX2

TILT-IN-SPACE WHEELCHAIR

Highly configurable and geared towards the rehabilitation environment, inclusive of pressure care cushioning surfaces, the Rehab RX2 is designed for users with more complex needs.



CONTOURED
BACKREST

CRASH TEST
APPROVED
AS/NZS 3696.19

1 YEAR WARRANTY
Components

2 YEAR WARRANTY
Frame



SWL
180kg



SWL
23.8st

Australia
1300 133 120 | aspirecare.com.au

United Kingdom
01622 541235 | aidacare.co.uk

Adjustable and Removable Headrest: Offers customisable support for the user's head and neck

Attendant-operated Drum Brakes and Push-to-lock Handbrakes: Provides both the user and attendant with safe and easy braking options

Removable, Height and Depth-adjustable Armrest : Allows for versatile arm support and easy transfers

Height & Depth-adjustable Backrest: Supports proper posture and comfort based on individual needs

Quick-release, Height and Depth-adjustable 24 inch Rear Wheels and 6 inch Castors: Enhances manoeuvrability and allows for easier transport and storage

Swing-away, Removable, Padded and Elevating Articulating Legrest: Ensures comfort and easy access for transfers or adjustments

Standard Anti-tippers: Adds an extra layer of safety to prevent tipping



View our Aspire Wheelchair Accessories Range




ADJUSTMENTS




Backrest Height



Seat & Backrest Depth



Elevating Legrest



Tilt-In-Space
26° Tilt | 30° Recline

REHAB RX2 WHEELCHAIR	Seat Width mm / in	Seat Depth mm / in	Front Seat Height mm / in	Rear Seat Height mm / in	Overall Width mm / in	Backrest Height mm / in	Armrest Height mm / in	Min Chair Weight kg / lbs	SWL kg / st
MWS449710	410 - 460/ 16 - 18	440 - 530/ 17.3 - 20.8	410 - 440/ 16.1 - 17.3	400 - 475/ 15.8 - 18.7	690/ 27.2	460 - 710/ 18 - 28	310 - 460/ 12.2 - 18.1	41/ 90.4	180/ 28.3
MWS449720	460 - 510/ 18 - 20	440 - 530/ 17.3 - 20.8	410 - 440/ 16.1 - 17.3	400 - 475/ 15.8 - 18.7	740/ 29.1	460 - 710/ 18 - 28	310 - 460/ 12.2 - 18.1	43/ 94.8	180/ 28.3
MWS449740	510 - 560/ 20 - 22	470 - 560/ 18.5 - 22.1	410 - 440/ 16.1 - 17.3	400 - 475/ 15.8 - 18.7	790/ 31.1	460 - 710/ 18 - 28	310 - 460/ 12.2 - 18.1	47.2/ 104.1	180/ 28.3

Note: Range of motion dimensions are accurate to a +/-10mm tolerance

ASFY1413 070425