



# Active Reach Tray Table

Style, comfort, and convenience. Always within reach.

A stylish, adjustable companion designed to bring everyday tasks closer and make life at home that little bit easier.



### Extendable Base

Flexible design to sit neatly with most sofa, chair, and recliner sizes.



### Adjustable Height

The smooth height adjustment system locks securely, allowing you to position the table exactly where you need it.



### Versatile Everyday Surface

The spacious, Carrara Marble-style laminate tabletop with removable cup holder provides an easy clean surface to maintain daily hygiene.



Carrara Marble Laminate

### Support to Stand and Sit

The integrated stand-assist handle offers a stable point of balance, helping you stand or sit with greater confidence.

### Simple, Space-saving Storage

When not in use, the tray table folds down to a slim silhouette, making it easy to slide under a bed, behind furniture, or into a cupboard without cluttering your space.

**1** YEAR WARRANTY

Code	Tray Height	Tray Size	Handle Height	Base Length*	Product Weight	Handle SWL	Tray SWL
BEA821815	660 - 810mm	400 x 450mm	870 - 1020mm	560 - 790mm	9.75kg	113kg	13.5kg

\*Suits chairs with distance between front and rear leg