



ADJUSTABLE DAY CHAIR - PRESSURE REDUCING

Comfort | Pressure Management | Infection Control



Pressure Redistribution

- Carrflex™ PU 4-Way Stretch Fabric provides superior chemical resistance whilst retaining flexibility and breathability.
- High quality Australian memory foam provides a comfortable and supportive surface, improving user comfort and pressure redistribution.
- Maximised Immersion and Envelopment which is fundamental for a support surface to effectively redistribute pressure improving blood flow which promotes healing.



FEATURES & BENEFITS

- **Contoured Seat and Backrest Padding** designed for pressure redistribution and user comfort
- **Head Cushion** height adjustable and removable head cushion for extra comfort
- **Sturdy design** for a generous weight capacity of 160kg, yet lightweight for easy manoeuvrability
- **Height Adjustable Legs** with easy push button adjustment for comfortable positioning and changing needs
- **Tool Free Adjustment Points** allow the chair to be customised to suit long term use or readily adapted to suit multiple users
- **Adjustable Seat Depth, Width and Height** allows customised configuration to a user's needs and circumstances
- **Unique Width Adjustable Frames** cater to the full spectrum of users shapes and sizes

1 YEAR WARRANTY

S.W.L
160KG

OVERALL WIDTH
(mm)



OVERALL DEPTH
(mm)



OVERALL HEIGHT
(mm)



SEAT DEPTH
(mm)



SEAT WIDTH
(mm)



SEAT HEIGHT
(mm)



BACKREST HEIGHT
(mm)



OVERALL WEIGHT
(kg)



**ADJUSTABLE
DAY CHAIR**
PRESSURE REDUCING
CHP208035

575 - 675

760 - 860

1070 - 1220

475 - 590

510 - 610

410 - 530

740 - 790

26

1300 133 120

www.aspirecare.com.au

