

## ALTITUDE RECLINER

The Aspire Altitude Recliner is an electric recliner chair designed to provide comfortable seating for long periods of time. It is a truly mobile and flexible seating solution, with 75mm castors facilitating smooth user transport plus a lithium-ion battery that can power the recliner all day when not within reach of power points.

## Its suitable applications include:

- Aged care facilities; •
- Rehabilitation facilities; •
- Hospitals; •
- Blood and treatment clinics; •
- Care in the home; and .
- Other independent in-home use •



Removable armrests

## Features & Benefits

- Use all day without plugging in the power with its • inbuilt battery (requires charging)
- Easy relocation of the chair between residents or facility locations with smooth rolling 75mm Tente castors (rear locking)
- Reduce shear and promote optimal skin microclimate with a multi-stretch vapour-permeable seating surface designed for users at risk of pressure injury



Can accomodate ROHO® cushions (or equivalent)



Optional Oncology Armrest CHP227715 (Pair)



- Medical grade vinyl on the outer chair surface for improved longevity and infection control
- Removable armrests facilitate safe lateral transfers, and optional accessories - oncology armrests and armrest table (see below)
- Enhanced comfort and immersive qualities with a • three-tiered backrest cushion



**Optional Armrest Table** CHP227720 (Single)



Transportation handle



\* Components: Fabric and mechanism.

ELECTRONICS

YEAR WARRANTY

COMPONENTS\*