



APOLLO SCOOTER

SMALL PULL APART

Stylish and Portable for Everyday Independence

The Aspire Apollo Pull Apart Small Scooter offers an ideal blend of safety, style and convenience for users seeking reliable mobility support at home and in the community. Featuring a vibrant red finish, users can also personalise their scooter to suit their style using seven additional interchangeable colour panels.





Personalised Style

Supplied in a vibrant red, with seven optional plug-and-play colour panel accessories to suit your style. *(Green option shown)*



Pull-apart Portability

Breaks down into five parts for safe lifting, transport, and compact storage.



LED Battery Status Indicator

Clear battery monitor helps you travel with confidence.



Electromagnetic Brakes

Provide reliable, controlled stopping for added user safety.



LED Headlight

Easily visible for nearby pedestrians and vehicles in low-light environments.



Seat Functionality


Features a swivel, height-adjustable seat to support easy side transfers and ride comfort.

Overall L x W x H	1090 x 550 x 910mm
Seat Size D x W	400 x 450 - 550mm
Seat to Footrest Height	370 - 420mm
Armrest Height	180mm
Min Turning Radius	1200mm
Ground Clearance	75mm
Max Climbing Angle	6°
Max Driving Speed	6km/h
Driving Range	12km
Wheels (puncture proof)	190mm (Front) 203mm (Back)
Battery	2 x 12V 12Ah (sealed lead acid)
Product Weight (excl. battery)	45kg
Heaviest Part	18.4kg
Battery Weight	2 x 4kg
Motor	270W, 2 Pole, Brushed
Controller	DY 50A
Brake System	Electromagnetic
Free Wheel Mode	Yes
Charger Size	24V 2A



SCT671500
Red Panels (standard)

Optional Colour Panel Accessories



SCA669360
White



SCA669365
Green



SCA669370
Blue



SCA669375
Pink



SCA669380
Purple



SCA669385
Yellow



SCA669390
Orange

*Scooter Frame/Motors/Axles/Brakes 2 Years | Battery 6 Months | Electronics 1 Year

2 YEAR WARRANTY*

SWL
136kg

1300 133 120 | aspirecare.com.au

ASFY176.10.300725