



ASSIST 3 WHEELCHAIR

Versatile self-propelling wheelchair designed to meet the diverse needs of Community, Hire, Hospital, and Aged Care settings. This reliable and comfortable wheelchair offers a range of accessories to support both individuals and caregivers.



2 YEARS WARRANTY
Frame

1 YEAR WARRANTY
Components

- **Robust & Reliable:** Built with a strong steel frame, removable swing-away leg rests, and quick-release self-propelling rear wheels, this transport solution is both practical and durable
- **Low Maintenance & Portable:** Solid rear tyres and a center cross brace design ensure low maintenance and portability, perfect for facilities, hire, or community use
- **Swing-Away and Removable Arm Rests:** Designed for effortless user transfers, enhancing transport and storage convenience
- **Versatile Performance:** Large, puncture-proof wheels handle varied terrains with ease
- **User-Friendly Rear Wheel Brake:** Features a large handle for easy access and operation



Folds for easy transport and storage



Heel loops on the foot rests



Swing away and removable foot rests



Quick release 24" solid rear tyres

ASSIST 3 WHEELCHAIR

	Seat Width mm / in	Overall Width mm	Seat Depth mm	Seat Height mm	Backrest Height mm	Armrest Height mm	Min. Chair Weight kg	SWL kg
MWS452200 350mm Wide	350 / 14	570	400	460 – 500	440	210	17.3	140
MWS452210 400mm Wide	400 / 16	620	400	460 – 500	440	210	17.5	140
MWS452220 450mm Wide	450 / 18	670	400	460 – 500	440	210	17.7	140
MWS452230 500mm Wide	500 / 20	720	400	460 – 500	440	210	17.9	140

*Range of motion dimensions are accurate to a +/-10mm tolerance

ASFY15.4 150824