BABY BASSINET

Aspire for... Safety, Nurture & Care















Reduce the effects of reflux using the Bassinet elevation position.

FEATURES & BENEFITS

- Hospital-grade castors: 100mm castors with lockable levers for enhanced durability, traction and safety.
- Heavy-duty Stainless Steel frame:
 Wipe clean surfaces adhere to stringent infection control procedures.
- Bassinet: Transparent and durable ABS bassinet provides cradled comfort for newborns.
- Ward label: Asset identification label allows for easily accessible equipment tracking between wards.
- **Storage:** 2x Polypropylene utility tubs for accessible storage.
- Mattress: Seamless foam mattress with multi-stretch medical grade Carrflex fabric and an in-built mesh extrusion for optimal breathability.

CODE	OVERALL LENGTH	OVERALL WIDTH	OVERALL HEIGHT (INC. BASSINET)	BASSINET LENGTH	BASSINET WIDTH	PRODUCT WEIGHT	SWL
PED490240	795 mm	450 mm	1050 mm	710 mm	360 mm	20 kg	30 kg