



Bedside Commode Range

Offering a safe and reliable bedside toileting solution for users with reduced mobility, helping minimise transfers to the bathroom while supporting comfort, confidence and independence in the home.



Features

- Height adjustable
- Raised padded armrests
- Easy-clean vinyl surfaces
- Non-slip rubber tips for stability
- Durable powder-coated steel frame
- Comfortable, supportive seating

Classic
Practical, lightweight design with removable commode bowl with lid

Deluxe
Fold-away padded seat with removable commode bowl with integrated carry handle

Champagne Vinyl

Slate Vinyl

1 YEAR WARRANTY

Code	Description	Overall Width	Overall Depth	Overall Height	Seat Width	Seat Depth	Seat Height	Backrest Height	SWL
BEC045005	Classic - Champagne Vinyl	560	600	900 - 1040	450	440	475 - 615	450	125kg
BEC045006	Classic - Slate Vinyl	560	600	900 - 1040	450	440	475 - 615	450	125kg
BEC045115	Deluxe - Champagne Vinyl	600	600	865 - 1015	460	450	410 - 565	440	160kg
BEC045116	Deluxe - Slate Vinyl	600	600	865 - 1015	460	450	410 - 565	440	160kg