

## VIDA WHEELCHAIR

The Aspire VIDA is a super lightweight transport wheelchair ideal for quick trips to the shops, appointments and catching up with friends. Ergonomic backrest and seat shape follow the contours of the user to maximise comfort and postural support. The 21 stone safe working load caters to a broad range of users.











## **FEATURES & BENEFITS**

- · Available as attendant propelled or self-propelled
- Lightweight frame design ensures easy lifting in and out of cars, ideal for carers
- Contoured and padded seating improves comfort and postural support
- · Attendant operated hand brakes with parking lock
- · Simple and intuitive folding for easy transportation
- · Removable leg rests and folding back for compact storage
- · Rear storage pocket for essentials
- · Anti-tippers available (WCA848275)
- · Lap belt available (WCA848300)

## Self - Propelled





## **Attendant Propelled**









	OVERALL LENGTH mm/in	OVERALL WIDTH mm/in	OVERALL HEIGHT mm/in	SEAT DEPTH mm/in	SEAT WIDTH mm/in	ANTERIOR SEAT HEIGHT mm/in	POSTERIOR SEAT HEIGHT mm/in	BACKREST HEIGHT mm/in	ARMREST HEIGHT mm/in	REAR WHEEL in	MIN. CHAIR WEIGHT kg/lbs	TOTAL CHAIR WEIGHT kg/lbs
	$\bigoplus_{\longleftrightarrow}$	$\leftrightarrow$	<b>1</b>		$\leftrightarrow$		100		1			
Attendant	995 /	655 /	930 /	420 /	450 /	495 /	445 /	430 /	200 /	13	8.8 /	9.6 /
Propelled	39	26	36.5	16.5	18	19.5	17.5	17	8		19.4	21.2
Self	1030 /	685 /	930 /	420 /	450 /	495 /	445 /	430 /	200 /	22	7.5 /	11.7 /
Propelled	40.5	27	36.5	16.5	18	19.5	17.5	17	8		16.6	25.8