



ELIO SCOOTER

SMALL PULL APART

Smart, Strong and Exceptionally Portable

Experience the perfect blend of strength, comfort and convenience with the Aspire Elio Small Scooter. Built for an active lifestyle, Elio's lightweight design, flexible battery options, and comfort features offer confident, independent mobility, perfect for shopping, social visits, or days out.





Light to Lift

At just 10.8kg, the heaviest part is easy for carers to lift and fits neatly in most car boots or storage spaces.



Ultra-lightweight Aluminium Frame

Easily disassembles into five compact parts for effortless transport and storage.



Additional Battery Options Available

Standard range battery up to 24km. Extended and airline-approved batteries also available (up to 35 km).



Small Turning Circle

Ideal for navigating through tight spaces such as hallways, lifts, or small shops without frequent repositioning or reversing.




Shock-absorbing Suspension

4-wheel independent suspension absorbs bumps from uneven surfaces for a smooth ride.



Intuitive Controls

Navigate with confidence using an intuitive dashboard featuring an LED battery display, speed controls, and a horn.



Smooth Transfers, Ultimate Comfort

Featuring a height-adjustable, swivelling padded seat, this scooter ensures easy side transfers every time.



Flexible Tiller Angle for Personalised Handling

Quickly adjusts to your ideal position, delivering a smooth and comfortable driving experience.

Overall L x W x H	1048 x 550 x 945mm
Seat Size D x W	400 x 450 - 550mm
Seat to Footrest Height	400 - 475mm
Armrest Height	180mm
Min Turning Radius	1200mm
Ground Clearance	60mm
Max Climbing Angle	6°
Max Driving Speed	6.4km/h
Driving Range	24km (35km optional)
Front/Rear Wheels (puncture proof)	190mm puncture proof
Battery	1 x 24V 20Ah (Lithium)
Product Weight (excl. battery)	28.6kg
Heaviest Part	10.8kg
Battery Weight	3.2kg
Motor	270W, 2 Pole, Brushed
Controller	PG 45A
Brake System	Electromagnetic
Free Wheel Mode	Yes
Suspension	Yes
Charger Size	24V 2A

*Scooter Frame/Motors/Axles/Brakes 2 Years | Battery 6 Months | Electronics 1 Year



SCT671700
Blue



SCT671710
White

Optional Batteries

SCA669400
Airline-approved battery (11.2km)

SCA669410
Extended range battery (35km)

2 YEAR WARRANTY*

