

FLEX SEAT WALKER

- · Super lightweight making it easy to manoeuvre and easy to transport
- Available in 6" (Champagne) and 8" (Matte Black) wheel variations
- Robust aluminium frame for improved durability and stability
- Height adjustable seat height to suit most users
- Soft padded seat for improved seating comfort
- · Adjustable and lockable brake handles





Champagne 6"



Black







FLEX
ADJUSTABLE
SEAT WALKER









570mm

22.4"

Overall



300mm

11.8"

Seat



360mm

14.1"

Seat



430-520mm

16.9-20.4"

Seat



710mm

27.9"

Backrest

Height



290mm

11.4"

Folded



6"

Castor



Product

Weight

6" Wheels Champagne Frame
WAF709670BK

WAF709660CH

6" Wheels
Champagne Frame


WAF709670E		
8" Wheels		
Black Frame		

730–830mm	730mm
28.7–32.6"	28.7"
780–970mm	820mm

30.7–38.1"







760mm 29.2"

8" 390mm 15.3" 8" 6.3kg 13.8lb

7.3kg

16lb