FOREARM WALKER

Aspire for... Mobility, Support & Comfort









Lockable handbrakes for superior safety and confidence in mobility.

FEATURES & BENEFITS

- Height and depth adjustable contoured forearm pads for optimal fit and postural assistance.
- Swivel front wheels for advanced manoeuvrability and control.
- Soft, padded, rubber moulded contoured forearm pads for enhanced comfort.
- Sturdy and lightweight aluminium frame provides stability and support for users up to 180kg.
- Foldable frame for trouble-free storage and transportation.

CODE	FOREARM PAD HEIGHT	OVERALL LENGTH	OVERALL WIDTH	FOLDED LENGTH	BRAKE HANDLE LENGTH	WHEEL SIZE	PRODUCT WEIGHT	SWL
WAF699200	1010 - 1200 mm	790 mm	640 mm	280 mm	450 - 520 mm	7.5"	16 kg	180 kg

