



Gabriel X5

Five Action Lift Recline Chair

The Aspire Gabriel X5 Lift Recline Chair has powered adjustments for tilt, headrest, lumbar, backrest, and leg elevation to support your individual comfort needs. Three level backrest massage allows you to relax with a massage level to suit you.



Dark
Grey

2 YEAR
WARRANTY



Petite



Large



Zero Gravity
Tilt-in-space



Back
Recline



Leg
Raise



Power
Lumbar



Power
Headrest



Seat
Heating



3 Function
Massage



Zero Gravity and Tilt-in-space: Elevates the legs above the head assisting in redistributing pressure and relieving discomfort



Integrated Leg Extension: Provides improved support under the lower legs for a wide range of users



USB-A & USB-C Charging: Stay connected with conveniently placed charging ports in both the armrest and handset for phones, tablets, and other devices



Large Button Hand Controller: Offers high visibility back-lit buttons that are easy to use. Hand control loops are located on both sides, to increase flexibility of setup for user



Three Level Backrest Massage: Gives users a choice of settings that feels most comfortable



Durable and Safe: Scratch-resistant materials and protective covers extend the life of the chair, and a built-in safety barrier keeps pets away from moving parts



Electric Headrest Adjustment: Offers customisable neck positioning at the press of a button

Adjustable Seat Heating: Integrated heating provides gentle warmth and improves comfort, perfect for year-round relaxation

Five Action Lift and Recline Functions: Includes independent back recline, lift, leg raise, lumbar support, and headrest adjustment

All-day Battery Backup: A lithium-ion battery ensures continued operation during power outages, supporting safety and independence

Code	Size	Seat Height	Seat Depth	Seat Width	Backrest Height	Overall Height	Overall Width
CHP227592	Petite	475mm	450mm	490mm	710mm	1100mm	850mm
CHP227597	Large	500mm	520mm	530mm	780mm	1170mm	900mm