











The Aspire Go Turner is a proven sit-to-sit transfer aid for quick and easy manual handling. This compact folding patient mover has a very narrow footprint, making it the perfect solution for transferring patients in tight spaces. The pull apart and folding design reduces storage space requirements and helps to transport the product if needed.



	Overall Length	Overall Height	Overall Width	Knee Pad Height from Footplate	Handle Height from Footplate	Footplate Height from Floor	Overall Unit Weight
							
Go Turner LSS390585	685mm 26.9"	1205mm 47.4"	570mm 22.4"	230-400mm 9 - 15.7"	875mm - 1145mm 34" - 45"	60mm 2.4"	18.9kg 41.6 lbs

Standing Slings - SWL 150kg			
LSS391132 - SML	LSS391133 - MED	LSS391134 - LRG	LSS391135 - XL

*Range of motion dimensions are accurate to a +/-10mm tolerance