ASFY69.1.14 181121

HOMECARE BATHROOM RANGE



Aspire for... Stability, Independence and Support



OVER TOILET AID

Adjustable height raised toilet seat allows for safe ingress and egress when toileting



SHOWER CHAIR

Provides a complete seating support for users with reduced mobility during showering



SHOWER STOOL

Allows for a comfortable bathing experience for users who require that little bit of extra support

SHARED FEATURES

- Durable high grip feet for stability and safety
- Adjustable legs to suit a wide range of users
- Raised armrests assist the user in transferring on and off





	SEAT DEPTH (mm)	SEAT WIDTH (mm)	SEAT HEIGHT (mm)	OVERALL DEPTH (mm)	OVERALL WIDTH (mm)	WIDTH BETWEEN ARMS (mm)	OVERALL HEIGHT (mm)
Over Toilet Aid BTT145025	410	340	405 - 560	480	520	475	595 - 750
Shower Chair BTS118025	340	410	440 - 570	340	525	465	750 - 880
Shower Stool BTS119035	340	410	440 - 570	340	525	465	590 - 720

Tested to AS/NZS3973.7:2009

