aspire HYBRID PATIENT SLINGS



FEATURES & BENEFITS

- Hybrid sling facilitates both Hygiene and General Purpose manual handling
- Intuitive shape allows ease of access for toileting and dressing
- \cdot Full coverage shoulder and head support
- Security belt and padded leg straps support the trunk and thighs
- Durable loop webbing designed not to tangle or 'cork-screw' during cleaning
- \cdot Embroidered unique serial numbers for traceability
- Service record labels to easily record and identify service history
- $\cdot\,$ Machine washable up to $85^\circ\,$



CODE	DESCRIPTION	EDGING COLOUR	SWL
LSS390830	Aspire Hybrid Patient Sling - Small	•	200kg
LSS390835	Aspire Hybrid Patient Sling - Medium	•	200kg
LSS390840	Aspire Hybrid Patient Sling - Large	•	300kg
LSS390845	Aspire Hybrid Patient Sling - Extra Large	•	300kg