



IDAHO DUAL ACTION 2 MOTOR

- Contemporary medium size
- Gentle tapered armrests suit the modern home environment
- Simultaneous adjustable back and leg support functions
- Integrated lift and recline functions for easy sit-to-stand
- Generous padding in backrest for improved lumbar support



Leg Recline



Back Recline



Optional Neck Pillow

CHP227171 | Charcoal Fabric

CHP227172 | Latte Vinyl



Charcoal Fabric



Latte Vinyl

Code	Colour	Seat Height	Seat Depth	Seat Width	Backrest Height	Overall Height	Overall Width	SWL
CHP227160	Charcoal	480mm	480mm	480mm	620mm	1100mm	800mm	150kg
CHP227165	Latte							