



MATISSE QUATTRO ACTION

4 MOTOR

The Matisse is a quattro action lift and recline chair, designed for therapeutic support, easy sit-to-stand assistance, and personalised adjustments—ideal for improved posture, circulation, and relaxation.



Santa Barbara Fog
Microsuede



Back Recline



Leg Recline



Power Lumbar



Power Headrest

2 YEAR WARRANTY

SWL
159kg

- **Available in Petite and Large Sizes** fit to the user's body shape, improving comfort, posture and safety when standing
- **Quattro Action Recliner** with back recline and lift, legrest, lumbar support and headrest adjustments
- **Integrated Lift and Recline Functions** to assist sitting and standing
- **Adjustable Headrest**, ideal for watching TV whilst reclined
- **Legrest Extends** to provide support under the lower legs, especially beneficial for taller people
- **Large Side Storage Pockets** ideal for books and newspapers
- **All-day Battery Backup**, a lithium-ion battery ensures continued operation during power outages
- **Backlit Hand Control** with integrated USB charging point

Code	Size	Seat Height	Seat Depth	Seat Width	Backrest Height	Overall Height	Overall Width	SWL
CHP227555	Petite	455mm	455mm	510mm	560mm	1015mm	915mm	136kg
CHP227560	Large	500mm	545mm	570mm	730mm	1065mm	965mm	159kg