



MONET DUAL ACTION 2 MOTOR ZERO GRAVITY

- Contemporary Large size
- Dual Action: Independent adjustable headrest, back and leg support functions
- Backlit hand control with integrated USB charging point.
- Backrest heating for improved climate control
- Zero gravity - Tilt in space pre-set positioning creates a feeling of weightlessness
- Comes with integrated lithium battery backup



Zero Gravity



Full Recline



Power Headrest

Colour



Canyon Ocean
Microsuede



Code	Seat Height	Seat Depth	Seat Width	Backrest Height	Overall Height	Overall Width	SWL
Canyon Ocean - CHP227540	495mm	535mm	560mm	570mm	1065mm	1065mm	159kg

ASBR26.5 200923