

# PORTO SINGLE ACTION

## 2 MOTOR ZERO GRAVITY SPACE SAVER



Zero Gravity



Space Saver



Full Recline

Colour



Fawn  
Soft Suede

- Designed to utilise limited space when fully reclined, without compromising on the overall seat size
- Head and armrest covers for improved longevity
- Zero gravity function redistributes sitting pressure, assists with oedema and relaxes muscles
- Designed for more petite users

Code	Seat Height	Seat Depth	Seat Width	Backrest Height	Overall Height	Overall Width	SWL
Fawn - CHP227200	490mm	500mm	490mm	690mm	1120mm	770mm	150kg