

REMBRANDT QUATTRO ACTION 5 MOTOR ZERO GRAVITY



- Contemporary Medium size
- Integrated wireless phone charger
- Seat and armrest heating for improved climate control
- Genuine shiatsu roller massage in backrest provides an exceptional seating experience
- Integrated footrest extension suits most user heights
- $\cdot \,$ Comes with integrated lithium battery backup

- Independent adjustable headrest back, leg and lumbar support functions
- Integrated lift and recline functions for easy sit-to-stand
- Memory pre-set positioning to capture the users most comfortable position with the press of a button
- Zero gravity tilt in space pre-set positioning creates a feeling of weightlessness

Col	our





Canyon Walnut Microsuede

Code	Seat Height	Seat Depth	Seat Width	Backrest Height	Overall Height	Overall Width	SWL
Canyon Walnut – CHP227550	495mm	535mm	560mm	600mm	1095mm	1015mm	159kg