

# RENOIR QUATTRO ACTION

## 4 MOTOR

- Available in three sizes: Small, Medium & Large
- Independent adjustable headrest, back, leg and lumbar support functions
- Ergonomic backrest design reduces pressure off the lower sacral area for an improved seating experience
- Large side storage pockets, ideal for books and newspapers
- Integrated footrest extension suits most user heights
- Integrated lift and recline functions for easy sit-to-stand
- Backlit hand control with integrated USB charging point
- Adjustable headrest, ideal for watching TV whilst reclined
- Features an all-day lithium ion battery backup to maximise safety and independence
- Four swivel castors allow the chair to be moved around with ease



Back Recline



Leg Recline



Power Lumbar



Power Headrest

### Colour



Dusk  
Microsuede

**2** YEAR WARRANTY

Code	Size	Seat Height	Seat Depth	Seat Width	Backrest Height	Overall Height	Overall Width	SWL
Dusk – CHP227575	Small	450mm	510mm	510mm	660mm	1010mm	750mm	130kg
Dusk – CHP227580	Medium	470mm	530mm	550mm	760mm	1060mm	800mm	130kg
Dusk – CHP227585	Large	470mm	530mm	590mm	760mm	1060mm	830mm	130kg