



SIGNATURE 2 SINGLE ACTION

1 MOTOR SPACE SAVER

The Signature 2 Space Saver is ideal for close wall positioning, as the chair moves away from the wall while reclining, making the best use of the user's living space.



Optional Neck Pillow

CHP197781 | Mink Warwick
CHP197783 | Lagoon Warwick
CHP197782 | Heritage



- **The Signature 2 is Available in Three Sizes Small, Medium, and Large**, each thoughtfully designed to deliver the perfect blend of comfort, postural support, and usability for a customised seating experience
- **The Space Saver Recline Chair** shifts forward while reclining, allowing for closer wall positioning
- **Single Motor** with a single button press for both lift and recline actions is intuitive and assists cognition
- **Adjustable Seat Depth** offers 25mm of depth adjustment, allowing taller users to achieve a better seated position
- **Integrated Footrest Extension** helps with fatigue and improves lower limb comfort
- **Easy To Move** four swivel castors included as standard allows the chair to be moved around with ease
- **All-day Battery Backup** with a lithium-ion battery ensures continued operation during power outages, supporting safety and independence
- **Backlit Hand Control** with integrated USB charging point

2 YEAR WARRANTY



Code	Colour	Size	Seat Height	Seat Depth	Seat Width	Backrest Height	Overall Height	Overall Width	SWL
CHP228100 CHP228110 CHP22812	Mink Warwick Lagoon Warwick Heritage	Small	480mm	460 – 485mm	480mm	580mm	980mm	740mm	145kg
CHP228130 CHP228140 CHP22815	Mink Warwick Lagoon Warwick Heritage	Medium	480mm	490 – 515mm	510mm	620mm	1050mm	780mm	145kg
CHP228160 CHP228170 CHP22818	Mink Warwick Lagoon Warwick Heritage	Large	500mm	510 – 535mm	540mm	650mm	1110mm	780mm	145kg