



# MAXI SHOWER STOOL



- **Super heavy duty** adjustable width shower stool
- **Soft padded seat and arms** provide postural support
- **Reinforced aluminium frame** accommodates bariatric users and is corrosion resistant
- **Height adjustable legs** with easy C-Clip adjustment points for changing needs

- **Raised armrests** for assistance in transferring on and off the stool
- **Durable non-slip feet** provide excellent grip and stability in wet conditions

**5** YEARS WARRANTY  
Corrosion

**1** YEAR WARRANTY  
Structure

**SWL**  
**310kg**

## MAXI SHOWER STOOL

	Overall Height mm	Overall Width mm	Overall Depth mm	Seat Height mm	Seat Width mm	Seat Depth mm	Between Arms mm
<b>BTS119030</b>	670 - 790	645 - 845	495	455 - 575	575 - 775	455	575 - 775

✓ DVA CONTRACTED

✓ AS / NZS Approved

✓ MAXI is registered in the ARTC, AUST R ARTC 306118

ASFY134 211024