

MONET DUAL ACTION 2 MOTOR ZERO GRAVITY

- Contemporary Large size
- Dual Action: Independent adjustable headrest, back and leg support functions
- Backlit hand control with integrated USB charging point.
- Seat and backrest heating for improved climate control
- Zero gravity Tilt in space pre-set positioning creates a feeling of weightlessness
- \cdot Comes with integrated lithium battery backup



Code	Seat Height	Seat Depth	Seat Width	Backrest Height	Overall Height	Overall Width	SWL
Canyon Ocean – CHP227540	495mm	535mm	560mm	570mm	1065mm	1065mm	159kg