



Monet X5

Five Action Lift Recline Chair

Enjoy the ultimate comfort and support of the Monet X5 Lift Recline Chair. A sophisticated, five action seating solution designed to adapt effortlessly to your posture, enhance relaxation throughout the day, and provide independence at home.



Zero Gravity
Tilt-in-space



Back
Recline



Leg
Raise



Power
Lumbar



Powered
Headrest



Seat & Backrest
Heating



Midnight
Blue

2 YEAR
WARRANTY





Integrated Lift and Recline Functions: Assists with sit-to-stand transfers, helping users feel confident and independent at home



Zero Gravity and Tilt-in-space: Elevates the legs above the head assisting in redistributing pressure and relieving discomfort



Powered Headrest: Ensures your neck is always supported comfortably



Integrated Leg Extension: Provides improved support under the lower legs for a wide range of users



Clinical and Care-friendly Controller: Includes a multi-function hand control for comprehensive positioning, or an optional simplified handset to assist users with reduced dexterity or cognitive needs (Code: CHP228910)



Durable and Safe: Scratch-resistant materials and protective covers extend the life of the chair, and a built-in safety barrier keeps pets away from moving parts



Integrated Seat and Backrest Heating: Helps to maintain a comfortable temperature, promoting relaxation throughout the seasons

All-day Battery Backup: The lithium-ion battery ensures continued operation during power outages



Five Action Recliner: Features independent back recline, legrest raise, lumbar support, an adjustable headrest and zero gravity tilt-in-space

Large Side Storage Pockets: Keep books, magazines or remotes tidy and within reach

Code	Size	Seat Height	Seat Depth	Seat Width	Armrest Height	Backrest Height	Overall Height	Overall Width	Overall Depth
CHP227542	Petite	455mm	455mm	490mm	200mm	700mm	1080mm	830mm	860mm
CHP227543	Large	500mm	550mm	550mm	200mm	750mm	1170mm	890mm	960mm