

MONTANA HD LIFT RECLINE CHAIR

- **True Bariatric Seating:** 300kg and 450kg options available.
- **Dual Action:** Independent control of legrest and lift/recline functions for infinite positioning, to achieve optimal posture for comfort and oedema management.
- **Configurable Backrests:** Optional alternative backrests are available on custom made chairs including laterally supportive Immersion Cushion and Waterfall designs. All cushions include rear zippers allowing foam/support adjustment for posture and gluteal shelf management.
- **Custom Sizing:** Capability to offer clients a custom-sized lift chair with the ability to specify the seat width, depth, height and armrest height. This option is available via the custom order form.
- **Durable Construction:** commercial grade foams and fabric for long duration sitting and reclining.



HUGE 160kg Legrest Weight Capacity

Designed for oedema management, the legrest caters to static loads of up to 160kg.



Optional Backrests

Available on custom made chairs.



Immersion Cushion



Waterfall

Standard backrest cushion shown in main large image

WARWICK



Barley Vinyl
Standard Colour



2 YEAR WARRANTY
Mechanism & Electronics

3 YEAR WARRANTY
Upholstery

LIMITED LIFETIME
Frame











HD300



HD450

MONTANA HD DUAL ACTION

SEAT HEIGHT (mm)	SEAT DEPTH (mm)	SEAT WIDTH (mm)	SEAT TO TOP OF BACK (mm)	OVERALL HEIGHT (mm)	OVERALL DEPTH (mm)	OVERALL WIDTH (mm)	DISTANCE FROM WALL (mm)
							
500	500	650	810	1270	1000	900	500
500	500	850	810	1270	1000	1100	500

CUSTOM HD300
CUSTOM HD450

Please refer to the Montana HD order form for custom sizes