

MONTANA HD DUAL ACTION

2 MOTOR



Optional Backrests

Available on custom made chairs.





Immersion Cushion

Waterfall



Full Recline



HUGE 160kg Legrest Weight Capacity*







Barley Vinyl



Custom Upholstery available

- · Dual Action: Back Recline and Leg Recline
- True Bariatric Seating: 254kg and 450kg options available
- **Custom Sizing:** Specify the seat width, depth, height and armrest height. This option is available via the custom order form
- Configurable Backrests: Optional alternative backrests are available on custom made chairs including laterall supportive Immersion Cushion and Waterfall designs
- HUGE 160kg Legrest Weight Capacity loads of up to 160kg*



Code	Seat Height	Seat Depth	Seat Width	Backrest Height	Overall Height	Overall Width	Legrest SWL*	SWL
Standard HD CHP227405	500mm	500mm	650mm	810mm	1270mm	900mm	160kg	254kg
Standard HD450 CHP227410	500mm	500mm	850mm	810mm	1270mm	1100mm	140kg	450kg
Custom HD Custom HD450	Please refer to the Montanna HD order form for custom sizes							