



# OSCAR SCOOTER

## MIDI

### Agile, Reliable, and Built for Everyday Use

The Oscar Scooter is designed for users who value practical mobility and confident handling. Its compact frame, responsive controls, and reliable power make it perfect for everyday trips, whether in local shops, or around your community.







### Extended Travel Range

Large lead-acid battery offers up to 39km on a single charge, reducing the need for frequent recharging.



### Plush, Supportive Seating

Features a deeply padded, adjustable captain's seat with suspension, and flexible tiller. Seat slides, swivels and reclines for perfect positioning.



### Smart Safety

Emergency stop button and automatic speed reduction while turning keeps you safe in unexpected situations.



### Comfortable Ride

Independent suspension absorbs bumps and reduces vibrations for comfort on uneven surfaces.



### Simple and Intuitive Controls

Straightforward analogue interface is easy for new and experienced users to operate.



### Compact and Agile Navigation

With a 1400mm turning circle, it easily manoeuvres in tight indoor spaces and busy community areas.

Overall L x W x H	1200 x 640 x 1300mm
Seat Size D x W	450 x 480 - 530mm
Seat to Deck Height	440 - 490mm
Seat to Ground Height	585 - 635mm
Armrest Height	170mm
Min Turning Radius	1400mm
Ground Clearance	65mm
Max Climbing Angle	10°
Max Driving Speed	10km/h
Driving Range	39km
Tyres	290mm   Pneumatic
Battery	2 x 12V 50Ah (Lead Acid)
Product Weight (excl. battery)	71kg
Battery Weight	30kg
Power   Motor	700W   4 Pole
Controller	RHINO 2 90A
Brake System	Electro-Mechanical
Free Wheel Mode	Yes
Charger Size	24V 5A

\*Scooter Frame/Motors/Axles/Brakes 2 Years | Battery 6 Months | Electronics 1 Year



SCT671211  
Cherry Red



SCT671212  
Steel Blue



ASFY184.3 190825