



OSLO SEAT WALKER

Sleek lines, elegant curves and eye-catching shades give the OSLO seat walker sophisticated poise. The aero-space grade aluminium frame is highly durable and ultra lightweight for easy manoeuvrability, folding and lifting, while tool-free adjustments can suit individual user needs.



- **Ultra-Lightweight Design:** Constructed with aerospace-grade aluminium for exceptional strength, the 6.25kg minimum weight makes it easy to lift and transport
- **Simple to Adjust:** Tool-free handle and backrest height adjustment accommodates users of varying stature, with push-button locking for secure customisation
- **Ergonomic Comfort:** Contoured, non-slip hand grips to reduce wrist strain and adjustable backrest for tailored support
- **All-Terrain Mobility:** Puncture-resistant PU wheels with swivel front castors and fixed rear wheels improve manoeuvrability
- **Portable & Compact:** One-click folding mechanism collapses the rollator into a slim profile (20cm width) for effortless car boot storage or travel
- **User-Centric Accessories:** Secure zippered storage bag with 5kg capacity and built-in cane holder

2 YEAR WARRANTY
Aluminium Frame

1 YEAR WARRANTY
Components

**SWL
136kg**

OSLO SEAT WALKER

	Overall Height mm	Overall Length mm	Overall Width mm	Seat Width mm	Seat Height mm	Seat Depth mm	Folded Width mm	Front Wheel mm	Rear Wheel mm	Min Weight kg	Total Weight kg	SWL kg
 WAF705320PU Purple	830 - 940	660	540	450	585	205	200	230	195	6.25	6.9	136
 WAF705320GR Green												