



Raised Toilet Seat

Adjustable with swing-back arms.

Regain confidence and independence in the bathroom.



Anti-slip Support Armrests

Textured armrests allow users to safely grip the armrests when sitting and standing from the Raised Toilet Seat for support.

Swing-back for Easy Transfers with Walking Aids

Swing-back mechanism allows armrest to be moved out of the way to assist transfers from wheelchairs, and walking aids.

Hygienic and Easy to Clean

Including a lid and components made from moulded plastic, it is easy to clean to minimise infection risks.

Universal Fit

Featuring depth and height adjustable mounting bolts, the Aspire Toilet Seat Raiser is versatile to suit almost any toilet.

Height Adjustable

Three height options to suit users of various heights or with mobility limitations, helping them regain confidence and remain independent at home.



Code	Overall Width	Overall Height	Seat Height	Distance between armrests	Weight
BTT164050	565mm	370mm	50 / 85 / 120mm	480mm	4.2kg