HS928 SCOOTER

MAXI

A robust MAXI class scooter that really will go the distance!

This heavy duty scooter is ideal for larger users and for traveling longer distances with a travel range up to 50km! The sporty style and numerous safety features ensures it is both comfortable and safe, and a real eye catcher.



Overall Length (mm)

(mm)

Width

Seat Width (mm)

Wheel Size (inch)

(b)

Speed

(mm)

Clearance

Turning Radius (mm) Travel Range (km) Battery Size (V)

ttery Product Weight ize incl. batteries V) (kg)



SCT671310 SILVER

1600 720 x 1280

560

15/15

10

110

1860

50


2 x 12V 100AH


178


* Scooter Frame 3 Years | Motors/Axles/Brakes 18 Months | Battery 6 Months



SCOOTER ACCESSORY RANGE




















Scooters	1	2	3	4	5	6	7	8	9
	Oxygen Tank Holder	Cane Holder	Crutch Holder	Sun Shade Canopy	Rear Shopping	Scooter Flag	Cupholder	Dust Cover	Rear Bucket
HS520									
HS589									
HS828	SCA6686000	SCA668610	SCA668620	SCA668650	SCA668630	SCA668640	_	_	_
HS898									
HS928									
Superlite									
Sonic	SCA668330	SCA668230	SCA668390	SCA668380	SCA668340	SCA668370	SCA668170	SCA668210	SCA668180
GT	3CA000330	3CA000230	3CA000330	3CA000300	3CA000340	3CA000370	3CA000170	3CA000210	3CA000100
Pathrider									