

SUPALITE SCOOTER

SMALL PULL-APART

Every day is a new adventure with the Aspire Boot Scooter Supalite.

An exceptionally stable and manoeuvrable 4-wheel travel scooter. Benefiting from a class-leading, ultralightweight design and the ability to be disassembled into lightweight pieces to travel anywhere, the Boot Scooter Supalite keeps up with your active lifestyle!



FEATURES & BENEFITS

- Disassembles into six lightweight pieces
- Compact, foldable seat with armrests
- · Easy-to-grip tiller adjustment knob
- Stable front frame-mounted seat post
- · Front convenience basket
- · Deluxe tiller with USB charger
- · Cup holder





Available in WHITE



Overall Length (mm) Overall Width (mm) Seat Front / Rear Width Wheel Size

Max Speed (km/h) Ground Tui

Turning Radius (mm) Travel Ba Range S (km)

Battery Size (V) Product Weight Heaviest incl. batteries component (kg) (kg)









8/8















SCT670616 WHITE

1022

480 x 880

440

6.5

31

1385

18.5

1 x 24V 3 12 AH (Lithium Ion)

34

10.7kg


^{*} Scooter Frame/Motors/Axles/Brakes 2 Years | Battery 6 Months



SCOOTER ACCESSORY RANGE




















Scooters	1	2	3	4	5	6	7	8	9
	Oxygen Tank Holder	Cane Holder	Crutch Holder	Sun Shade Canopy	Rear Shopping	Scooter Flag	Cupholder	Dust Cover	Rear Bucket
HS520									
HS589									
HS828	SCA668600	SCA668610	SCA668620	SCA668650	SCA668630	SCA668640	_	_	_
HS898									
HS928									
Superlite									
Sonic	SCA668330	SCA668230	SCA668390	SCA668380	SCA668340	SCA668370	SCA668170	SCA668210	SCA668180
GT	3CA000330	3CA000230	3CA000390	3CA000300	3CA000340	3CA000370	3CA000170	3CA000210	3CA000100
Pathrider									