




# SHOWER STOOL WITH ARMS



BTS119000



BTS119025

- Clip on seat or padded back attachment for easy assembly and cleaning
- Easy push button height adjustable legs
- All frames individually treated for corrosion resistance
- Seat with anterior tilt to facilitate standing from sitting and provides greater comfort for stiff or painful hips
- Raised armrests for assistance in transferring on and off the stool

- Durable non-slip feet provide excellent grip and stability in wet conditions and a stainless steel washer enhances durability
- Additional height increments on one leg to accommodate uneven surfaces
- Tested to AS/NZS3973.7:2009

**5 YEARS WARRANTY**  
Corrosion

**1 YEAR WARRANTY**  
Structure



## SHOWER STOOL WITH ARMS

|                                      | Overall Height<br>mm | Overall Width<br>mm | Overall Depth<br>mm | Seat Depth<br>mm | Seat Width<br>mm | Seat Height<br>mm | Between Arms<br>mm | Product Weight<br>kg |
|--------------------------------------|----------------------|---------------------|---------------------|------------------|------------------|-------------------|--------------------|----------------------|
| STEEL<br>BTS119000                   | 670 - 855            | 545                 | 400                 | 310              | 575 - 775        | 505 - 655         | 495                | 3.8                  |
| ALUMINIUM<br>BTS119005               | 670 - 855            | 545                 | 400                 | 310              | 575 - 775        | 505 - 655         | 495                | 3.1                  |
| WITH PADDED SEAT<br>BTS119020        | 670 - 855            | 545                 | 400                 | 310              | 575 - 775        | 505 - 655         | 495                | 4.2                  |
| WITH PADDED SEAT & BACK<br>BTS119025 | 805 - 955            | 545                 | 400                 | 310              | 575 - 775        | 505 - 655         | 495                | 6.8                  |

