



SIERRA SCOOTER

MEDIUM

Smooth, Stable, and Ready for Everyday Adventures

This scooter combines a classic design with smooth, quiet operation and the power, travel range, and stability required for outdoor use over varying terrains.





Stable and Secure on Every Trip

4-wheel configuration with front and rear suspension provides a steady ride over varied terrain.



Extended Daily Travel

High-capacity lead-acid batteries deliver up to 40km of range per charge, reducing the need for frequent recharging.



Designed for Easy Manoeuvring

A tight 1565mm turning circle makes navigating narrow paths, shop aisles, and small community spaces simple.



Plush, Supportive Seating

Features a deeply padded, adjustable captain's seat and flexible tiller, seat slides, swivels and reclines for perfect positioning.



Simple and Intuitive Controls

Straightforward analogue interface is easy for new and experienced users to operate.




Smart Safety

Emergency stop button and automatic speed reduction while turning keeps you safe in unexpected situations.

Overall L x W x H	1290 x 625 x 1150mm
Seat Size D x W	450 x 480 - 530mm
Seat to Deck Height	440 - 490mm
Seat to Ground Height	585 - 635mm
Armrest Height	215mm
Min Turning Radius	1565mm
Ground Clearance	87mm
Max Climbing Angle	12°
Max Driving Speed	10km/h
Driving Range	40km
Tyres	290mm Pneumatic
Battery	2 x 12V 50Ah (Lead Acid)
Product Weight (excl. battery)	67kg
Battery Weight	30kg
Power (Nominal Peak)	700W 1400W
Motor	4 Pole
Controller	RHINO 2 120A
Brake System	Electro-Mechanical
Free Wheel Mode	Yes
Charger Size	24V 5A

*Scooter Frame/Motors/Axles/Brakes 2 Years | Battery 6 Months | Electronics 1 Year



SCT671261
Arctic Mist



SCT671271
Desert Sand



2 YEAR WARRANTY*

SWL
136kg

ASFY189.9 151225