

# SIGNATURE 2 SPACE SAVER

## Space Saver

The Signature 2 Space Saver is ideal for closer wall positioning as the chair tracks away from the wall as it reclines, conserving the best use of your living space. The simple handset with single button press for both lift and recline functions is intuitive and easy to use, ideal for those with reduced cognitive ability. Tailored with premium Warwick Global Suede and Vinyl upholstery and high grade durable foams to provide years of postural support, comfort and safety.



*Pictured in Mink*

## Features & Benefits

- **Simple handset** with single button press for both lift and recline actions is intuitive and assists cognition.
- **Space Saver** - Chair shifts forward during recline action to allow closer wall positioning.
- **Integrated lift and recline functions** for easy sit-to-stand and maximum comfort.
- **Integrated lithium-ion battery backup** ensures chair can be used in the event of a power outage or if moved between rooms during the day.
- **Depth adjustable seat** caters to taller users.
- **4 x swivel castors** included as standard, allows the chair to be moved around with ease.
- **Backlit hand control** with integrated USB charging point.
- **Headrest and armrest protectors** included as standard

### WARWICK



**Mink**  
Microsuede



**Lagoon**  
Microsuede



**Heritage**  
Vinyl

**2 YEAR WARRANTY**


**S.W.L 145KG**

SIGNATURE 2 SPACE SAVER	SEAT HEIGHT (mm)	SEAT DEPTH (mm)	SEAT WIDTH (mm)	SEAT TO TOP OF BACK (mm)	OVERALL HEIGHT (mm)	OVERALL DEPTH (mm)	OVERALL WIDTH (mm)	DISTANCE FROM WALL (mm)
<b>SMALL</b>								
Mink CHP228100	480	460-485	480	580	980	690	740	160
Lagoon CHP228110								
Heritage Vinyl CHP228120								
<b>MEDIUM</b>								
Mink CHP228130	480	490-515	510	620	1050	715	780	160
Lagoon CHP228140								
Heritage Vinyl CHP228150								
<b>LARGE</b>								
Mink CHP228160	500	510-535	540	650	1100	740	780	160
Lagoon CHP228170								
Heritage Vinyl CHP228180								

# SIGNATURE 2 DUAL ACTION

## Dual Action

The Signature 2 Dual Action offers more positioning, for increased comfort and support. With independent control of back and legrest functionality, the user has infinite positioning options. The smart legrest extends in length as it rises, supporting more of the lower legs and enhancing comfort particularly for taller users. Tailored with premium Warwick Global suede and vinyl upholstery and high grade durable foams to provide years of postural support, comfort and safety.



*Pictured in Lagoon*

## Features & Benefits









- **Dual Action** - Independent control of leg and recline for infinite positioning.
- **Extending legrest** supports more of the lower legs for enhanced comfort.
- **Depth adjustable seat** caters to taller users.
- **Integrated lithium-ion battery backup** ensures chair can be used in the event of a power outage or if moved between rooms during the day.
- **4 x swivel castors** included as standard allows the chair to be moved around with ease.
- **Backlit hand control** with integrated USB charging point.
- **Headrest and armrest protectors** included as standard.

### WARWICK



**2 YEAR WARRANTY**

**S.W.L 155KG**

SIGNATURE 2 DUAL ACTION	SEAT HEIGHT (mm)	SEAT DEPTH (mm)	SEAT WIDTH (mm)	SEAT TO TOP OF BACK (mm)	OVERALL HEIGHT (mm)	OVERALL DEPTH (mm)	OVERALL WIDTH (mm)	DISTANCE FROM WALL (mm)
								
<b>SMALL</b>								
Mink CHP228200	480	460-485	480	580	980	710	740	480
Lagoon CHP228210								
Heritage Vinyl CHP228220								
<b>MEDIUM</b>								
Mink CHP228230	480	490-515	510	620	1050	735	780	480
Lagoon CHP228240								
Heritage Vinyl CHP228250								
<b>LARGE</b>								
Mink CHP228260	500	510-535	540	650	1100	760	780	480
Lagoon CHP228270								
Heritage Vinyl CHP228280								