



TRANSIT 3 WHEELCHAIR

The Aspire TRANSIT 3 wheelchair is a durable transport chair designed for use by an attendant or carer. It is suitable for community, hospital, or aged care settings.

Available in durable Dacron fabric, or in hospital-grade vinyl with fully welded seams and concealed mounts to meet strict infection control standards.



1 YEAR WARRANTY
Components

2 YEAR WARRANTY
Frame

Transit 3

SWL
140kg

Transit 3 WARD Vinyl HD

SWL
180kg

INFECTION CONTROL WARD VINYL VERSION

- ✓ Durable vinyl that is easy to wipe clean
- ✓ Fully welded seams to inhibit bacterial growth
- ✓ Concealed mounting of upholstery to prevent fluid ingress

Folding Steel Frame: Durable and compact for easy storage and transport

Removable Legrest & Footrest: Swing-away for easy transfers, height adjustable for comfort

Removable Armrests: Padded and easily removable for smoother transfers

12" Quick Release Wheels & Adjustable Castors: Ensures smooth movement, with quick release for easy storage

Push-to-lock Brakes: Locks the chair securely for stability during stops

Padded Seat & Backrest: Comfortable for short-term use

Fire Retardant: Safe and durable Dacron or Vinyl upholstery

Lightweight: Easy to handle and transport

Anti-tippers: Included as standard (WCA848270)



View our Aspire Wheelchair Accessories Range

TRANSIT 3 WHEELCHAIR

	Seat Width mm / in	Overall Width mm	Seat Depth mm	Anterior Seat Height mm	Posterior Seat Height mm	Backrest Height mm	Armrest Height mm	Min. Chair Weight kg	SWL kg
400mm Dacron MWS452300	400 / 16	540	400	460 - 500	450	440	210	16.4	140
450mm Dacron MWS452310	450 / 18	590	400	460 - 500	450	440	210	16.6	140
500mm Dacron MWS452320	500 / 20	640	400	460 - 500	450	440	210	16.8	140
450mm WARD Vinyl MWS449445	450 / 18	610	400	460 - 500	450	440	210	16.5	140
500mm WARD Vinyl MWS449450	500 / 20	660	400	460 - 500	450	440	210	16.8	140
570mm HD WARD Vinyl MWS449480	570 / 22	730	400	460 - 500	450	440	210	17.9	180

*Range of motion dimensions are accurate to a +/-10mm tolerance

ASFY118.5 070425