



Velaro & Velaro Plus

Walking Sticks


The Velaro & Velaro Plus Walking Stick features a lightweight aluminium frame, 12-level height adjustment, and an ergonomic contoured handle for enhanced comfort and support.



Gloss Black



Gloss Purple



Matte Green



Matte Orange






Velaro Walking Stick




T-handle



Single Tip



Velaro Plus Walking Stick



Swan Neck



Multi-point Tip

12-level Height Adjustment



Loosen

Turn collar left



Adjust

Press pin and slide until it clicks into place



Secure

Turn collar right until snug



Gloss Black

Velaro WAS748700


Velaro Plus WAS748800



Gloss Purple

Velaro WAS748710

Velaro Plus WAS748810



Matte Green

Velaro WAS748720

Velaro Plus WAS748820



Matte Orange

Velaro WAS748730

Velaro Plus WAS748830

Description	Overall Height	Product Weight	SWL
Velaro	685 – 964mm	0.4kg	136kg
Velaro Plus	730 – 1010mm	0.6kg	136kg