



VIDA X WHEELCHAIR

The Aspire VIDA X is an ultra-modern, super-lightweight transport wheelchair designed for mobility convenience, comfort and style.

The wheelchair's ergonomic backrest and seat shape meticulously follow the user's contours, providing maximum comfort and optimal postural support.


With 22 inch sports mag wheels and eye-catching frame & upholstery detail, the VIDA X is a striking modern transport wheelchair solution.



UPHOLSTERY
Black
FRAME
Matte Black
MWS449867



UPHOLSTERY
Black/Purple
FRAME
Shiny Purple
MWS449868



LIGHTWEIGHT, COMFORTABLE AND STYLISH

- Self-propelled
- Lightweight aluminium body for easy handling and portability
- Enhanced comfort seating via padded seat and backrest with subtle upholstery detail
- Sporty 22 inch mag wheels
- Soft padded armrests
- Swing-away legrests are height adjustable for enhanced flexibility and comfort
- Storage pouch neatly stores essential items
- Choose the matte black frame with striking red upholstery stitching or opt for shiny purple frame with funky purple stitching
- Folds down for easy transport and convenient storage
- Anti-tippers available (WCA848275)
- Lap belt available (WCA848300)



Folds for easy transport and storage



Sporty rear wheels



Modern frame colours and upholstery detailing



Storage pouch

ASPIRE VIDA X WHEELCHAIR	Overall Length	Overall Width	Overall Height	Seat Depth	Seat Width	Anterior Seat Height	Posterior Seat Height	Backrest Height	Armrest Height	Minimum Chair Weight	Total Chair Weight	Rear Wheel
	mm/ in	mm/ in	mm/ in	mm/ in	mm/ in	mm/ in	mm/ in	mm/ in	mm/ in	kg/ lbs	kg/ lbs	in
	1040 / 41	690 / 27	950 / 26	420 / 16.5	450 / 18	500 / 19.6	460 / 18	430 / 17	200 / 8	7.5 / 16.6	12.2 / 26.7	22

*Range of motion dimensions are accurate to a +/-10mm tolerance

ASFY104.T1 220724