



VIKING WALKING TUTOR

The Viking Walking Tutor is designed for hospital and clinical use, providing stable, comfortable, and fully adjustable support for upright walking. Its strong aluminium frame, padded arm supports, and easy-to-use hydraulic height adjustment make it ideal for rehabilitation and assisted mobility in care settings.



Optional Accessory
WAA648340 - Bag
WAA648350 - Holder



Optional IV Pole
WCA871190

- **Hydraulic Height Adjustment:** Lever-operated system makes it easy to raise or lower the arm supports smoothly and safely, helping users of different heights
- **Adjustable Padded Arm Supports:** Soft armrests offer comfortable forearm support, for users with limited grip strength. The width can be adjusted to suit individual body shapes and provide added stability
- **Width-adjustable Base:** Improves manoeuvrability in tight spaces and around clinical equipment
- **Angle-adjustable Soft Handgrips:** Supports a natural wrist position, making it easier and more comfortable for users to maintain control
- **Optional Accessories:** Available to enhance functionality

1 YEAR WARRANTY

**SWL
200kg**

VIKING
WALKING TUTOR
WAF699130

Overall Height	Overall Length	Overall Width	Width Between Forearm Supports	Wheel Size	Product Weight	Bag SWL
900 - 1325mm	850mm	730 - 890mm	520 - 750mm	100mm	18.3kg	5kg