



# Vogue Light Seat Walker

**Lightweight Aluminium Body** Improves manoeuvrability and portability

This modern European designed seat walker has a clean and crisp frame design and is available in a selection of eye-catching colours and prints to really stand-out while keeping you on the move.

Less Than  
**6kg**



**Handbrakes Lock**  
For user safety and stability



**Compact Turning Circle**  
For improved mobility in tight spaces



Standard



**Leopard**  
WAF705310LE



**Flower**  
WAF705310FL



**Zebra**  
WAF705310ZB



**Red**  
WAF705310RE



**Champagne**  
WAF705310CH

Tall




**Lime**  
WAF705315GR



**Ergonomic Handles** With personalised height adjustment to promote comfortable and correct posture

**Deep Seat and Wide Backrest** For improved comfort and support  
**Bonus** Walking cane holder



**Adjustable Handles**  
Allows for personalised positioning, helping to maintain proper posture



**Simple Folding Mechanism**  
Locks for easy transportation and storage



**Removable Storage**  
Zip bag secures up to 5kg of essentials



**Accessories**  
With a cane holder included as standard, choose from our range of optional accessories

|          | Overall Length<br>mm | Overall Width<br>mm | Overall Height<br>mm | Seat Width<br>mm | Seat Height<br>mm | Folded Width<br>mm | Wheel Size<br>mm | Minimum Functional Weight<br>kg | Typical Functional Weight**<br>kg |
|----------|----------------------|---------------------|----------------------|------------------|-------------------|--------------------|------------------|---------------------------------|-----------------------------------|
| Standard | 655                  | 620                 | 780 - 920            | 465              | 540               | 220                | 200              | 5.4                             | 6.45                              |
| Tall     | 655                  | 620                 | 805 - 1025           | 460              | 630               | 220                | 203              | 5.86                            | 6.6                               |

Range of motion dimensions are accurate to a +/-10mm tolerance | \*\*Includes accessories, bag, and back strap.