




Watson

Adjustable Pressure Reducing Chair

Designed to deliver exceptional comfort, configurability, and safety across a wide range of care environments. Built with a robust steel frame and industry-leading adjustability, it is an ideal solution for users requiring pressure seating that can be adjusted to suit diverse needs.



Pressure Injury Prevention

Supportive memory foam and four-way PU stretch fabric assists with immersion and pressure redistribution during extended use



High-SWL Steel Frame: Robust steel frame achieves a Safe Working Load (SWL) of 200kg, making it ideal for heavy-duty and long-term care environments

Optimised Postural Support: Contoured cushioning combined with adjustable head, seat, and back support promotes proper alignment and sustained comfort

Exceptional Adjustability: Fully adjustable seat width, depth, and height to accommodate varied body shapes, sizes, and postural needs

Tool-free Adjustment Points: Ensure quick configuration for multiple users

Integrated Spring-pin Safety System: Prevents over-extension and maintains secure, controlled adjustments

Multi-layered Commercial-grade Foam: Enhances longevity and sitting comfort

Mobility-ready Rear Castors: Helps transport the chair between rooms



Code	Seat Width mm	Seat Depth mm	Seat Height mm	Backrest Height mm	Armrest Height mm	Overall Depth mm	Overall Width mm	Overall Height mm
CHP208040	500 - 600	400 - 570	460 - 570	800	200	810 - 990	690 - 790	1160 - 1260