

WEIGHT BEARING PLATFORM SHOWER COMMUNE



Attendant Propelled Model



Self Propelled Model



View our Aspire Shower
Commode Accessory Flyer

- **Weight Bearing Platform** easy-slide weight bearing platform footrest slides completely away for transfers
- **Angle Adjustable Backrest** from 8° to 18° for improved seating comfort
- **Ultra Low Seat Height** of 590mm which can be reduced to 490mm when customised
- **Toilet Clearance** improved over toilet clearance and 8° backrest angle ensures patient comfort
- **Integrated Access Trough** improves hygiene whilst maintaining toilet bowl clearance height of 465mm
- **Open Armrests Style** ensures wider sitting surface and smaller base footprint for improved accessibility
- **Safety Arms** prevent falls and swing back out of the way for transfers
- **Brake Access** 4 x **TENTE** total locking castors provide brake access from all angles

Commode is registered in the ARTG AUST R ARTG 286695

5 YEAR WARRANTY
Corrosion

1 YEAR WARRANTY
Components



Code	Description	Height	Width	Seat Width	Backrest Height	Armrest Height	Anterior Seat Height	Posterior Seat Height
BTC066000	460mm Attendant Propelled	995	535	460	305	225	590	560
BTC066200	460mm Self Propelled	995	665	460	305	225	590	560
BTC066050	530mm Attendant Propelled	995	605	530	305	225	590	560
BTC066055	530mm Self Propelled	995	735	530	305	225	590	560

All measurements are in millimetres unless otherwise specified