

BED RAIL PROTECTORS

Made of protective soft foam, the Lifecomfort Bed Rail Protectors are designed to prevent patient injury from contact with bed rails.

Easy to apply and remove to any bed rail, making them applicable to most bed types.

Covers are durable and easy to clean, ensuring infection control standards can be adhered to.



- Durable and Easy To Clean Covers
 - Easy To Apply and Remove
 - Does Not Impede Bed Rail Function



0



BED RAIL PROTECTORS



Patient Protection Protects patients from injury by preventing limbs from protruding through bed rails as well as any contact



Velcro Fastening Straps Easy to apply and remove



Functionality Does not impede routine bed rail functionality



Easy Clean PU Fabric Covers are easily wipe cleaned to assist with infection control management

SPECIFICATIONS

Size	Small	Large
Unit	Pair	
Description	Bed Rail Protector	
Code	BEA010775	BEA010776
Width (folded)	335mm	490mm
Width	670mm	980mm
Length	1260mm	1580mm
Thickness	25mm	25mm



Sold in Pairs Sold in pairs and available in 2 sizes for optimal fit

CARE

Washable up to 95°c with water and mild detergent

Do Not use solvents or non-approved bleaches

Do Not iron

Do Not twist or wring out